

prietest™ Clinical Chemistry Reagents
GAMMA-GT

In vitro diagnostic test kit, for professional use only


INTENDED USE : Quantitative in vitro determination of activity of gamma-glutamyltransferase (Gamma-GT) in serum on photometric systems.

ORDERING INFORMATION	Pack Size	Cat No.
	1 X 10 ml	GGT 01 10
	1 X 25 ml	GGT 01 25
	2 X 50 ml	GGT 02 50
	4 X 50 ml	GGT 04 50

CLINICAL SIGNIFICANCE : Gamma-glutamyltransferase (GGT) is a membrane-localized peptidase mainly present in kidneys, pancreas, liver and prostate. This enzyme plays a significant role in glutathione metabolism and takes part in the transport of amino acids into the cells. The rise of GGT activity, often isolated (earlier and longer increase compared to other enzymes), is one of the most sensitive indicators of an affection of the liver or bile ducts. The strongest increases are observed in cases of intra hepatic or posthepatic biliary obstructions (reaching levels from 5 to 30 times normal), primary or secondary neoplasms, acute or chronic pancreatitis, and other pancreatic malignancies (especially those associated with hepatobiliary obstruction). More moderate elevations are observed during infectious hepatitis, cirrhosis and hepatic steatosis. Alcohol in chronic ingestion, some drugs like phenobarbital and phenytoine or contraceptives can also increase GGT rate in the serum.

METHOD :
Modified IFCC method. Enzymatic. Kinetic, UV.

PRINCIPLE :
Optimized kinetic determination of γ -glutamyltransferase (γ -GT) activity :



The increase of absorbance at 405 nm due to the formation of 5-amino-2-nitrobenzoic acid is proportional to γ -GT activity.

REAGENTS :

COMPONENTS AND CONCENTRATIONS:

N.B: Concentrations are those in the final test mixture.

R1 :

Tris Base	100 mmol/l
Glycyl Glycine	100 mmol/l

R2 :

L-Gamma-glutamyl-3-carboxy-4-nitroanilide	4 mmol/l
Preservative & Stabilizer	

STORAGE INSTRUCTIONS AND REAGENT STABILITY :

The reagents are stable up to the end of the indicated date of expiry on the vial label, if stored at 2 to 8°C, protected from light and contamination is avoided. Do not freeze the reagent! Discard the Reagent if found turbid or in case the Absorbance of Working Reagent exceeds 0.8 AU at 405 nm against distilled water.

WARNINGS AND PRECAUTIONS :

Take the necessary precautions for the use of laboratory reagents.

WASTE MANAGEMENT :

Please refer to local regulation requirements.

REAGENT PREPARATION :

- Two Reagent procedure : The reagents are ready-to-use.
- One Reagent Procedure : Mix four volumes of reagent R1 with one volume of reagent R2. Stability of working reagent solution : Three weeks at 2 to 8°C

MATERIAL REQUIRED BUT NOT PROVIDED :

NaCl solution 9 g/l, General laboratory equipment, Analyser / Photometer, pipettes etc.

SPECIMEN :

Serum free from hemolysis.
 Stability in serum : Sera are stable 7 days at 2 to 8°C, Two months at -20°C.
 Discard contaminated specimens.

ASSAY PROCEDURE 1: Two Reagent procedure

Application sheets for automated systems are available on request.

Wavelength	: Hg 405 nm
Optical path	: 1 cm
Temperature	: 37°C
Mode	: Kinetic

Bring all the contents of the kit to Room Temperature prior to use.
 Read rate of change of absorbance of sample against distilled water or air.
 Label the test tube as sample, control and pipette into respective test tube the reagent, sample, control sample as per the table given below :

	Sample / Control
Reagent R1	800 μ l
Reagent R2	200 μ l

Mix and pre warm the working reagent at 37°C for 2 minutes then add

Sample / Control	50 μ l

Mix and after a 60 seconds incubation at 37°C measure the change of absorbance per minute (ΔA /minute) during 180 seconds.

ASSAY PROCEDURE 2: One Reagent procedure

Label the test tube as sample, control and pipette into respective test tube the reagent, sample, control sample as per the table given below :

Prewarm working reagent at 37°C for two minutes prior to addition of sample.

	Sample / Control
Working Reagent	1000 μ l
Sample / Control	50 μ l

Mix and after a 60 seconds incubation at 37°C measure the change of absorbance per minute (ΔA /minute) during 180 seconds.

CALCULATION :

At 405 nm with one reagent procedure and two reagent procedure for 1 cm path light cuvette

Activity of Sample (U/L) = (ΔA /Min) X 2211

TEMPERATURE CONVERSION FACTORS :

To correct result to other temperatures multiply by factor shown in table

Assay Temperature	Conversion Factor to		
	25°C	30°C	37°C
25°C	1.00	1.37	1.79
30°C	0.73	1.00	1.30
37°C	0.56	0.77	1.00

CALIBRATION :

For the calibration of automated photometric systems use of the commercially available calibrator is recommended.

QUALITY CONTROL :

To ensure adequate quality, use of the commercially available control sera is recommended.

PERFORMANCE CHARACTERISTICS :

MEASURING RANGE :

The test has been developed to determine Gamma GT activity within a measuring range from 5 to 1000 U/L. When values exceed higher limit of the range, such samples should be diluted 1 + 1 with NaCl solution (9 g/l) and the result multiplied by 2.

SPECIFICITY/INTERFERENCE :

No interference was observed by Ascorbic Acid up to 20 mg/dl, Bilirubin up to 20 mg/dl, Glucose up to 500 mg/dl and lipemia up to 600 mg/dl Triglycerides. A list of drugs and other interfering substances with Gamma-GT determination has been reported by Young et al.

SENSITIVITY/ LIMIT OF DETECTION :

The lower limit of detection is 5 U/L.

PRECISION :

Intra-assay precision n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	25	0.2	0.80
Sample 2	60	0.5	0.83
Sample 3	336	3.0	0.89

Inter-assay precision n = 20	Mean [U/L]	SD [U/L]	CV [%]
Sample 1	28	0.2	0.71
Sample 2	63	0.5	0.79
Sample 3	344	3.0	0.87

METHOD COMPARISON :

A comparison between Robonik Prietest Gamma GT (y) and a commercially available test (x) using 20 samples gave following results:

Linear Regression : $y = 0.9956x - 3.0$ U/L

Correlation Coefficient : $r = 0.9990$

REFERENCE RANGE :

Men : 0 to 55 U/L **Women** : 0 to 38 U/L

It is recommended that each laboratory should assign its own reference range.

LITERATURE :

- Scherwin, J.E, *Liver function*. Clinical Chemistry: Theory, Analysis, Correlation, 4th Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 492 and appendix.
- Henderson, A.R., Moss, D.W., *Enzymes*, Tietz Fundamentals of Clinical Chemistry, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), 2001, 352.
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, 1995

INSTRUMENT APPLICATION prietest TOUCH	PARAMETERS FOR INSTRUMENT SETTING	
Name : GAMMAGT , Mod : KIN	TEST NAME	GAMMA GT
Pri.: 405 , Sec.: 0	Reaction	Kinetic
Temp: 37C , KF : 1.000	Reaction Slope	Increasing
Vol : 500ul , Unit : U/L	Wavelength 1	405 nm
Lag : 60 , Read : 180	Temperature	37°C
Blk : N, QC : Y, Norm : Y	Zero Setting	Distilled Water
Std : N , Factor : 2211	Factor	2211
Normal HI = 55	Units	U/L
Normal LO = 0	Sample Volume	50 μ l
QCNH : *	Reagent Volume	1000 μ l
QCNL : *	Lag Time	60 Seconds
QCABH = *	Read Time	180 Seconds
QCABL = *	Reference Range	0 to 55
Init. OD : = 0.8 H	Reagent Linearity	1000
Max Delta/Min : = 0.452	Max Delta/Min	0.452
Rgt. Linearity : 1000	Initial OD	< 0.8
NOTE :		
* Indicates user definable parameter.		
NA Implies Not Applicable		

prietest is the Trade Mark of ROBONIK (INDIA) PVT.LTD., for Clinical Chemistry Reagents.

prietest TOUCH is the Trade Mark of ROBONIK (INDIA) PVT. LTD., for Biochemistry Analyser.



An ISO 9001 : 2008 Certified Company
 An ISO 13485 : 2012 Certified Company

Manufactured and Marketed by:

ROBONIK (INDIA) PVT. LTD.,
 A-374, TTC Industrial Area, Mahape,
 Navi Mumbai - 400 710, INDIA.

Tel. No.: + 91 (22) 67829700

Fax. No.: + 91 (22) 67829701

Orders : sales.reagents@robonikindia.com

Queries : feedback@robonikindia.com

Website : www.robonik.in